

Table 1

Test Case ID	Test Scenario	Test Case	Pre-condition	Test Steps	Test Data	Expected Results
TC01	Test audio personalisation	Execute audio playback with audio personalisation(with guest account)	None	1. Open Mimi App on Android TV	Year of birth: 1982 Sound personalisation: active Intensity: 75%	Play audio with personalisation with guest account; successfully pause and play; disable and enable
				2. Navigate to settings and select		
				3. Select “Guest Mode”		
				4. Put in the date stated in Test Data		
				5. Go back to play mode		
				6. Activate personalisation and click play		
				7. Disable personalisation		
				8. Enable personalisation		
TC02	Test audio personalisation	Execute audio playback without audio personalisation(exit guest account)	Run after TC01	1. Click Settings		Personalisation is deactivated; The activation button is no more light green
				2. Click “Exit guest mode”		
				3. Confirm		
				4. Click play		
TC03	Switch to another video	Play another video using remote, while current video is playing		1. Open Mimi App on Android TV		Navigate through videos; Selected video plays with our error
				2. Click arrow down on remote		
				3. Select video by clicking right button 2 times		
				4. Click select button		
TC04	Switch to another video	Play another video using remote(no arrow down)	Run after TC03 and pause video	1. Click left arrow		First selected video is playing; Second selected video is playing
				2. Click select button		
				3. Click right arrow		
				4. Click select		
TC05	Login(no guest mode)	Login via QR code	1.Preinstall Mimi hearing app on mobile and create an account	1. Click arrow let and navigate to settings	Hearing Test account in mobile	The audio personalisation is activated The “activate personalisation” button is bright/ active
			2. Open Mimi app in TV	2. Select settings		
				3. Click “Continue”		
				4. Click “Login”		
				5. Scan QR Code with		
				6. Success page is shown		
				7. Go Bach to video		

Test Case ID	Actual Results	Status(Pass/Fail)
TC01		Pass
TC02	Personalisation is deactivated; The activation button is no more light green	Pass
TC03	Navigate through videos; Selected video plays with our error	Pass
TC04	The focus is changing between the control on the Botton tray; You need to click the bottom arrow(down) to be able to switch between videos	Fail
TC05	The video is still playing; Personalisation is active	Pass